Oracle Health Sciences Empirica Signal Upgrade Guide





Oracle Health Sciences Empirica Signal Upgrade Guide, Release 9.1.0.1

F37234-01

Copyright © 2020, Oracle and/or its affiliates.

This software and related documentation are provided under a license agreement containing restrictions on use and disclosure and are protected by intellectual property laws. Except as expressly permitted in your license agreement or allowed by law, you may not use, copy, reproduce, translate, broadcast, modify, license, transmit, distribute, exhibit, perform, publish, or display any part, in any form, or by any means. Reverse engineering, disassembly, or decompilation of this software, unless required by law for interoperability, is prohibited.

The information contained herein is subject to change without notice and is not warranted to be error-free. If you find any errors, please report them to us in writing.

If this is software or related documentation that is delivered to the U.S. Government or anyone licensing it on behalf of the U.S. Government, then the following notice is applicable:

U.S. GOVERNMENT END USERS: Oracle programs (including any operating system, integrated software, any programs embedded, installed or activated on delivered hardware, and modifications of such programs) and Oracle computer documentation or other Oracle data delivered to or accessed by U.S. Government end users are "commercial computer software" or "commercial computer software documentation" pursuant to the applicable Federal Acquisition Regulation and agency-specific supplemental regulations. As such, the use, reproduction, duplication, release, display, disclosure, modification, preparation of derivative works, and/or adaptation of i) Oracle programs (including any operating system, integrated software, any programs embedded, installed or activated on delivered hardware, and modifications of such programs), ii) Oracle computer documentation and/or iii) other Oracle data, is subject to the rights and limitations specified in the license contained in the applicable contract. The terms governing the U.S. Government's use of Oracle cloud services are defined by the applicable contract for such services. No other rights are granted to the U.S. Government.

This software or hardware is developed for general use in a variety of information management applications. It is not developed or intended for use in any inherently dangerous applications, including applications that may create a risk of personal injury. If you use this software or hardware in dangerous applications, then you shall be responsible to take all appropriate fail-safe, backup, redundancy, and other measures to ensure its safe use. Oracle Corporation and its affiliates disclaim any liability for any damages caused by use of this software or hardware in dangerous applications.

Oracle and Java are registered trademarks of Oracle and/or its affiliates. Other names may be trademarks of their respective owners.

Intel and Intel Inside are trademarks or registered trademarks of Intel Corporation. All SPARC trademarks are used under license and are trademarks or registered trademarks of SPARC International, Inc. AMD, Epyc, and the AMD logo are trademarks or registered trademarks of Advanced Micro Devices. UNIX is a registered trademark of The Open Group.

This software or hardware and documentation may provide access to or information about content, products, and services from third parties. Oracle Corporation and its affiliates are not responsible for and expressly disclaim all warranties of any kind with respect to third-party content, products, and services unless otherwise set forth in an applicable agreement between you and Oracle. Oracle Corporation and its affiliates will not be responsible for any loss, costs, or damages incurred due to your access to or use of third-party content, products, or services, except as set forth in an applicable agreement between you and Oracle.

Contents

Preface

Documentation accessibility	
Related resources	
Access to Oracle Support	
Additional copyright information	
Apply the 9.1.0.1 Patch to Oracle Empirica Signa	ıl
Apply the 9.1.0.1 Patch to Oracle Empirica Signa Before you start	
	ı l 1 1
Before you start	1
Before you start Prepare to upgrade	1



Preface

This preface contains the following sections:

- Documentation accessibility
- Related resources
- Access to Oracle Support
- Additional copyright information

Documentation accessibility

For information about Oracle's commitment to accessibility, visit the Oracle Accessibility Program website at http://www.oracle.com/pls/topic/lookup?ctx=acc&id=docacc.

Related resources

All documentation and other supporting materials are available on the Oracle Help Center.

Access to Oracle Support

Oracle customers that have purchased support have access to electronic support through Support Cloud.

Contact our Oracle Customer Support Services team by logging requests in one of the following locations:

- English interface of Oracle Health Sciences Customer Support Portal (https:// hsgbu.custhelp.com/)
- Japanese interface of Oracle Health Sciences Customer Support Portal (https:// hsgbu-jp.custhelp.com/)

You can also call our 24x7 help desk. For information, visit http://www.oracle.com/us/support/contact/health-sciences-cloud-support/index.html or visit http://www.oracle.com/pls/topic/lookup?ctx=acc&id=trs if you are hearing impaired.

Additional copyright information

This documentation may include references to materials, offerings, or products that were previously offered by Phase Forward Inc. Certain materials, offerings, services, or products may no longer be offered or provided. Oracle and its affiliates cannot be held responsible for any such references should they appear in the text provided.



1

Apply the 9.1.0.1 Patch to Oracle Empirica Signal

Before you start

You can apply the Oracle Empirica Signal 9.1.0.1 patch to an instance of Oracle Empirica Signal 9.1.0.0.360.

Prepare to upgrade

To upgrade, you must confirm that you have the appropriate credentials.

Unpack the Oracle Empirica Signal 9.1.0.1 files
 On the application server running Oracle Empirica Signal 9.1.0.0.360, start a command shell.

Deploy Oracle Empirica Signal

After you have unpacked the 9.1.0.1 files on the application server running Oracle Empirica Signal, you must deploy the Signal application in the Oracle WebLogic console.

Confirm the upgrade
 Confirm that the upgrade succeeded and notify users.

Before you start

You can apply the Oracle Empirica Signal 9.1.0.1 patch to an instance of Oracle Empirica Signal 9.1.0.0.360.

For information about installing Oracle Empirica Signal 9.1, see the Oracle Empirica Signal 9.1 *Installation and Upgrade Instructions*.

Prepare to upgrade

To upgrade, you must confirm that you have the appropriate credentials.

For this procedure, you must have credentials for the following accounts:

- An Oracle Empirica Signal administrator.
- A non-privileged user who can log in to the application server.
- An Oracle WebLogic Server administrator.
- 1. Using the administrator account, log in to the Oracle Empirica Signal application.
- To confirm that the version number of Oracle Empirica Signal is 9.1.0.0.360, click Settings, then About.
- 3. Click Back, and then View Currently Logged in Users.

If sessions in addition to your own are running, instruct the users to finish their sessions.

4. Click Exit.

Unpack the Oracle Empirica Signal 9.1.0.1 files

On the application server running Oracle Empirica Signal 9.1.0.0.360, start a command shell.

- 1. As the non-privileged user, on the application server running Oracle Empirica Signal 9.1, start a command shell (for example, Putty).
- 2. Create a temporary directory <Signal9101> (for example, /u01/stage/ Signal9101), and unpack ESI_9.1.0.1.zip into it. For example:

```
$ unzip ESI_9.1.0.1.zip -d /u01/stage/Signal9101
```

- 3. In the temporary directory, open the file create_Signal-9_1_0_1_361.sh and verify that directory /u01/stage/Signal_Install, specified in the file, exists. Otherwise, update the INSTALL_BASE_DIR, INSTALL_DIR, and SIGNAL_DIR variables with the actual path.
- 4. From the temporary directory, give the shell script execute permissions and run it:

```
$ chmod +x create_Signal-9_1_0_1_361.sh
$ ./create_Signal-9_1_0_1_361.sh
```

Verify that the script prints Completed and that there are no errors.

- 5. Verify that there is a backup copy of <INSTALL_DIR> (e.g., Signal_Install.9_1_0_0_360_bkp). Verify the <INSTALL_DIR>/Signal/WEB-INF/lib directory includes the following files:
 - Signal-9 1 0 1 361.jar
 - Datamine-1_5_1_1_361.jar
 - EmpiricaCore-1_5_1_1_361.jar
 - Topics-9_1_0_1_361.jar

Deploy Oracle Empirica Signal

After you have unpacked the 9.1.0.1 files on the application server running Oracle Empirica Signal, you must deploy the Signal application in the Oracle WebLogic console.

- 1. Log in to the WebLogic Server Administration Console.
- 2. On the left, click Lock & Edit.
- 3. From the panel on the left, select the **Deployments** link.
- 4. In the table on the right, select the **Signal** checkbox.
- Click Update, then Next, and then Finish.
- 6. On the left, click Activate Changes.
- 7. On the left panel, expand **Environment** and select **Servers**.
- 8. On the right, click the **Control** tab, then select the **SignalServer** checkbox.
- Click Shutdown/Force Shutdown now.
- 10. On the left, above the table, click the **Refresh** icon.



Wait until SignalServer status is SHUTDOWN.

- 11. Select the SignalServer checkbox.
- 12. Select the Start button.
- 13. Click Yes.
- **14.** If Topics Service is enabled, repeat steps 8 to 13 for **TopicsService** and **TopicsServer**.
- 15. Click the Refresh icon.

Wait until the statuses of SignalServer and TopicsServer (if enabled) are RUNNING.

Confirm the upgrade

Confirm that the upgrade succeeded and notify users.

- 1. Using the administrator account, log in to Oracle Empirica Signal.
- 2. From the left navigation pane, select **Settings**, then **About**.
- 3. Verify that the Oracle Empirica Signal version is 9.1.0.1.361.
- 4. Click Exit.
- 5. Notify users that the Oracle Empirica Signal instance is available.

